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

# Preface

Please read the user manual carefully before using the product.

1. This manual contains operating instructions, assembly methods, and simple solutions to malfunctions.
2. This manual applies to our model: 9004 G-MOVE C2.
3. This manual contains maintenance and self-check methods for the wheelchair. Keep it in a suitable place.
4. Provide this manual as a reference when other people are going to use this wheelchair.
5. The annotations and illustrations in this manual may differ slightly from the actual parts due to quality improvements or design changes. The actual product prevails.
6. Please contact your dealer if you have any questions or concerns.
7. Improper use of any vehicle can result in injury. Unsafe driving can cause harm to yourself and others.
8. The electric wheelchair is designed to comfortably transport people who have difficulty walking or are unable to walk.
9. This electric wheelchair is designed to transport one person only.

## SYMBOLS USED IN THIS MANUAL

	Product reference
	Medical device
	CE certificate
	Manufacturer
	Distributor
	Multiple use in the same patient
	Date of manufacture
	Batch number
	Instructions for use
	Caution/Warning
	Importer

## COMPONENTS OF THE CHAIR

This product complies with the standards set out in EU and UK regulations. Optional accessories and equipment are available at an additional cost.

# Index

<b>1. User information</b>	65
1.1 This instruction manual	65
1.2 For further information	65
1.3 Packaging	65
<b>2. Product presentation</b>	66
2.1 Regulations	66
2.2 Intended use	66
2.3 Product performance	67
2.4 Product information	68
2.5 Technical characteristics	69
2.6 Correct use	69
2.7 Contraindications	69
<b>3. Safety instructions</b>	70
3.1 Main safety features	70
3.2 General precautions when driving	70
3.3 Safety: Temperature	71
3.4 Safety: Moving parts	71
3.5 Safety: Lifting the wheelchair	71
3.6 Safety: Risk of suffocation	71
3.7 Other warnings	72
3.8 Service life	72
3.9 Electromagnetic compatibility	72
<b>4. Use and operation</b>	73
4.1 Assembly and adjustment of the wheelchair	73
4.2 Checks before using the wheelchair	75
4.3 Handling the wheelchair	76
4.4 Curves	76
4.5 Brakes and emergency stop	76
4.6 Obstacles and kerbs	77
4.7 Manual mode of the wheelchair	77
<b>5. Battery and charger</b>	78
5.1 Charger requirements	78
5.2 Using the charger	79
5.2.1 Type 1 (direct charging of the battery)	79
5.2.2 Type 2 (charging via the controller)	79
5.2.3 Type 3 (indirect battery charging)	79
5.3 Battery use and maintenance	80
<b>6. System diagnostics</b>	80
<b>7. Controller</b>	82
7.1 Control panel	82
7.2 Using the controller	82
<b>8. Maintenance</b>	84
8.1 Battery cleaning and maintenance	84
8.2 Wheelchair faults and checking	85
8.3 Maintenance checks	85
<b>9. Warranty</b>	86
9.1 Warranty issues	86
<b>10. Warning</b>	87
<b>11. Other</b>	87
<b>Appendix to the report: Manufacturer's declaration on EUT</b>	88

# 1. Information for the user

Thank you very much for purchasing a García 1880 wheelchair.

Our products are constantly being improved, so they may change without prior notice. However, any changes to the information provided by existing users must be clearly communicated if safety is seriously affected.

All technical specifications are approximate and may be subject to change.

The estimated service life of this product is 5 years. Please do NOT use or install third-party components in the wheelchair unless they have been officially approved.

## 1.1 This instruction manual

This manual will assist you in safely using and maintaining your wheelchair.

Do not use your wheelchair until you have read and understood the entire contents of this manual (including all booklets).

If any of the user manuals are not included with your wheelchair, please contact your distributor immediately.

## 1.2 For further information

Please contact your dealer if you have any questions about the use, maintenance, or safety of your wheelchair.

## 1.3 Packaging

In order to be eligible for a product refund or to obtain approval for air transport of your wheelchair, if applicable, you must keep the original packaging in which the product was supplied. You will need to repack your wheelchair in the original box to comply with air transport certification.

**García 1880** - Sucesores de Manuel García S.L  
Ctra Segorbe Castellnovo Km 0  
12400 Segorbe - Castellón  
Spain  
[www.garcia1880.com](http://www.garcia1880.com)

## 2. Product presentation

### 2.1 Regulations

This product is a Class I medical device and complies with Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices.

This product must not be disposed of in the waste bin. There are specific collection points where you can deposit them. For the disposal of the packaging and the product, strictly comply with the legal regulations of your community.

### 2.2 Intended use

- Technical aid to facilitate the mobility of persons with disabilities or reduced mobility, enabling them to remain seated and be transported over very short distances.
- Wheelchairs are designed exclusively for the personal use of a user who cannot walk or has reduced mobility, to be self-propelled or moved by a third party (assistant pushing the chair) indoors and outdoors.
- The product must be used under the specified conditions and for the intended purposes; otherwise, the warranty will be void.
- Driving a wheelchair requires cognitive, physical and visual skills. The user must be able to assess and correct the results of actions taken while operating the wheelchair.
- The wheelchair cannot carry more than one person at a time.
- The maximum weight limit (which includes the weight of the user plus the weight of any other accessories installed on the wheelchair) is indicated in this manual and on the product label.
- The user must familiarise themselves with the contents of this manual before operating the wheelchair. In addition, before beginning to operate the wheelchair, the user must have received complete instructions from a qualified authorised dealer. The first few sessions with the wheelchair should be conducted under the supervision of a qualified technician/trainer/dealer/caregiver.

## 2.3 Product performance

This series of electric wheelchairs is powered by a lithium-ion battery and a direct current motor. Users control the direction and adjust the speed using a joystick. It is suitable for low speeds, good road conditions and small slopes.

## 2.4 Product information

### PERFORMANCE PARAMETERS

Key technical data	9004		
Maximum load capacity	100 Kg	Engine	24V 220W brushless 2PCS
Total depth with footrest	930 mm		
Overall width	600 mm	Controller	35A
Total height	950 mm	Battery	/
Folded size	600*370*700mm		/
Packaging size	665*390*800mm		20Ah(3.5kg)
Net weight	19.4 kg		/
Gross weight	24.8 kg		/
Seat width	440 mm		/
Seat depth	430 mm		/
Seat height	460 mm		Charger
Backrest angle	5°	Maximum speed	≤ 5.5 km/h
Backrest height	530 mm	Horizontal braking distance on roads	≤ 1.0 m
Distance by which the headrest protrudes from the backrest	/	Safe braking on maximum slope	≤ 1.6 m ( 3°)
Headrest height	/	Retention performance on slopes	6°
Distance from footrest to seat	380 mm	Static stability	≥ 6°
Height of footrest from floor	95 mm	Dynamic stability	≥ 3°
The angle of the footrest relative to the seat surface	12°	Obstacle clearance height	≤ 25 mm
Distance between armrests	440 mm	Width of cross trench	100 mm
Front wheel diameter	180 mm/ 7 inch	Climbing capacity	<= 6°
Rear wheel diameter	220 mm/ 8,6 inch	Minimum turning radius	≤900 mm
Service life	5 years	Theoretical driving distance	≥ 20 km

## 2.5 Technical characteristics

- **Material:** Carbon fibre
- **Motor:** 2 x 24V 220W Silent, brushless
- **Battery:** Lithium 20 Ah (3.5 kg)
- **Maximum speed:** 5.5 km/h
- **Maximum distance travelled:** 20 km
- **Climbing capacity:** 6°
- **Turning radius:** 90 cm
- **Maximum kerb height:** 2.5 cm
- Solid wheels
- Adjustable anti-tip wheels included
- The backrest is not foldable
- Joystick configurable for right- and left-handed users
- **Electromagnetic brakes:** If the joystick is not used, the chair brakes automatically
- **Safety belt:** Yes
- **Horn:** Yes
- **LED light:** Yes
- **Remote control:** No
- **Joystick for companion:** No
- **Seat and backrest covers:** Removable and machine washable
- **Storage bag:** No

## 2.6 Correct usage

- The maximum weight limit (which includes the weight of the user plus the weight of any other accessories installed on the wheelchair) is indicated in this manual and on the product label.
- Do not use on motorways.
- Wheelchairs must always be used for their intended purpose. They must not be used for any other purpose. For example, they must not be used as toys, reclining chairs, to transport more than one person, etc.
- This electric wheelchair must never be carried by the armrests, as this could damage the structure of the chair. Check that the armrests are correctly positioned.
- The footrest must be at least 65 mm above the ground to prevent possible falls.
- Do not make any modifications (mechanical alterations) that could compromise your safety when using the wheelchair.
- Do not exceed the maximum weight.
- Do not drive under the influence of alcohol or when extremely fatigued.

## 2.7 Contraindications

The use of a wheelchair is prohibited if any of the following conditions prevent its safe operation:

- Visual disturbances or other perceptual disturbances.
- Altered consciousness or alertness.
- Uncontrolled involuntary movements or uncorrectable imbalance.
- Unsteadiness when sitting.

Please note that operating a wheelchair requires sufficient cognitive, physical and visual abilities. The user must be able to assess the effects of their actions while operating the wheelchair and, if necessary, correct them. These abilities and the safe use of additional attachments are not covered by García 1880's assessments. The manufacturer accepts no liability for any damage resulting from this.

### **WARNING!**

- If you are under the influence of any medication that could affect your ability to drive, you should not operate the wheelchair.
- In order to operate the wheelchair safely, the user must have adequate visual ability.
- The wheelchair cannot accommodate more than one person at a time.
- Do not allow children to ride in the wheelchair without supervision.

## 3. Safety instructions

### 3.1 Main safety features

- Classified according to type of protection against electric shock: **Internal power supply.**
- Classified according to type of protection against electric shock: **Type B application.**
- Classified according to degree of protection against liquid ingress: **IPX3.**
- Classified according to safety of use in a flammable anaesthetic mixture with air or gas mixed with oxygen or nitrous oxide, flammable anaesthetic gases: **Non-AP/APG type.**
- Classified according to mode of operation: **Continuous operation.**
- Nominal voltage: **DC 24V.**
- Does not provide protection against defibrillator shocks in that section.
- No part emits or receives electronic signals.
- Equipment not permanently installed.

### 3.2 General precautions when driving

- Keep your hands on the armrests to control the wheelchair.
- Practise in parks or other open, safe places until you can handle the wheelchair with skill.
- Practise driving thoroughly in a safe place to master the principles of moving forward, stopping and turning.
- If you are driving in the dark, wear light-coloured or reflective clothing to be more visible.
- To avoid falls and dangerous situations, it is recommended that you start using your new electric wheelchair on a flat surface with good visibility.
- Before driving on the road, have a companion or carer with you and check that it is safe.
- Strictly follow the traffic rules for pedestrians; do not consider yourself a vehicle driver.
- Drive on the pavement and zebra crossings. Do not drive in the vehicle lane.
- Drive smoothly, avoid driving in a zigzag pattern or turning sharply.
- If the wheelchair gets wet, dry it immediately.
- Avoid contact with sea water: sea water is caustic and can damage the wheelchair.
- Avoid contact with sand: sand can get into the moving parts of the wheelchair and cause excessive wear.
- Do not use the wheelchair if the temperature is below -25 °C or above +50 °C.
- Do not push or pull objects while using the wheelchair.
- Do not drive through puddles of water.

#### AVOID DRIVING IN THESE CONDITIONS

- Driving in adverse weather conditions, such as rainy days, thick fog, strong winds, snow, etc.
- Driving on poor road surfaces, such as muddy roads, trails, sandy roads, gravel, etc.
- Driving over unfenced ditches, ponds, etc.
- Crossing railway tracks.
- When you have to cross railway tracks, stop at the crossing to confirm that it is safe and make sure that the wheels do not get stuck on the tracks.
- The electric wheelchair is for personal use only; do not transport people or goods, and do not use it for traction purposes.

## PRECAUTIONS FOR GOING UP AND DOWN SLOPES

- Avoid driving in the following locations: steep hills, sloping areas, high steps, channels, etc.
- Avoid driving on steep slopes; the slope angle should be less than 6 degrees (Figure 1). Handle the controller carefully when driving on slopes.
- Reduce speed when going downhill.
- Avoid driving sideways on steep slopes.


Figure 1

Do not switch the wheelchair to manual mode when going up or down slopes. If the wheelchair malfunctions at a crossroads, ask someone passing by for help. Switch the wheelchair to manual mode and push it out of the way, or move to a safe place immediately.

## PRECAUTIONS FOR CARERS

- Caregivers should confirm that the user's feet are in the correct position on the footrests and ensure that clothing does not get caught in the wheels.
- Caregivers should push the wheelchair forward and maintain safety on steep or long slopes.

## 3.3 Safety: Temperature

- Avoid physical contact with the wheelchair motors at all times. The motors are continuously in motion during use and can reach high temperatures. After use, the motors cool down slowly. Physical contact can cause burns. The motors need at least 30 minutes to cool down.
- When not in use, ensure that the wheelchair is not exposed to direct sunlight for extended periods of time. Certain parts of the wheelchair, such as the seat, backrest, and armrests, may become hot if exposed to direct sunlight for too long. This could cause burns or allergic skin reactions.

## 3.4 Safety: Moving parts

All wheelchairs have moving and rotating parts. Contact with moving parts can cause serious physical injury or damage to the wheelchair. You must avoid contact with the moving parts of the wheelchair.

- Wheels (small and central).
- Footrest.
- Armrests.
- Latch for folding the backrest.
- Foldable control panel support.

## 3.5 Safety: Wheelchair lift

Do not lift the seat system by holding onto the removable parts, as this could damage it or cause injury to the user.

## 3.6 Safety: Risk of suffocation

This mobility aid uses small parts that, in certain circumstances, may pose a choking hazard to young children.

### 3.7 Other notices

#### REPAIR AND RECONDITIONING

If the wheelchair needs to be repaired or reconditioned, contact the seller or technical service department. Do not modify it yourself. This could cause an accident or malfunction of the wheelchair.

#### KEEP IT DRY

Do not place the wheelchair in damp locations. If the wheelchair gets wet, dry it immediately.

Do not modify the assembly or materials of this wheelchair without confirmation from our company.

To avoid causing imbalances, do not add weight arbitrarily.

When someone is seated in the wheelchair or the clutch is not in manual mode, do not use other vehicles to pull or push the wheelchair.

### 3.8 Service life

The shelf life of this product is 5 years from the date of manufacture. Use the product within the specified period. Do not use it more than one year after the expiry date to avoid accidents.

**Date of manufacture:** see label.

### 3.9 Electromagnetic compatibility

The wheelchair should be kept away from strong magnetic fields and large inductive electrical equipment, such as radio stations, television stations, underground radio stations, and radio stations that transmit mobile phones. Pay attention to whether there are sources of electromagnetic interference nearby and, if so, move away from them to avoid electromagnetic interference. The electric wheelchair must avoid electromagnetic interference.

Note:

- Electric wheelchairs must comply with the electromagnetic compatibility requirements of standard YY0505.
- The user must install and use the electric wheelchair based on the electromagnetic compatibility information provided.
- Portable and mobile RF communication devices may affect the performance of the electric wheelchair. Therefore, avoid strong electromagnetic interference, such as proximity to mobile phones or microwaves.
- Please refer to the appendix for the manufacturer's notice and declaration.
- The battery cut-off voltage is 23 V.
- The electric wheelchair belongs to class D of the GB/T 18029.21-2012 standard, a wheelchair with electronic differential steering and manual brakes.

The electric wheelchair must not be placed or used in conjunction with other devices. If necessary, check that the electric wheelchair can operate normally in these circumstances.

## 4. Use and operation

### 4.1 Adjusting the wheelchair

#### UNFOLDING THE WHEELCHAIR

Hold the backrest with one hand and the seat with the other, and open it with force (image 1). Fully unfold the wheelchair, then secure the lock located under the backrest (image 2). Before using it, ensure that the lock is securely fastened. Otherwise, it could fold while driving and cause a hazard.

Image 1

Image 2

#### FOLDING WHEELCHAIR

Loosen the fastener (image 3), then hold the backrest with one hand and pull the seat with the other hand, folding it down (image 4).

Image 3

Image 4

#### INSTALL THE DRIVER

Insert the controller into the armrest tube. Please note that you must keep the controller in a horizontal position (image 5). Then tighten the screw (image 6).

Image 5

Image 6

## BATTERY INSTALLATION AND REPLACEMENT

### Remove the battery

Press and hold the orange button and remove them (image 7).

Image 7

## FOOTREST PLATFORM

The platform can be dismantled (image 8).

Image 8

## ADJUST THE ANTI-TIP WHEELS

Turn the red nut (image 9).

Image 9

## FOLDABLE ARMREST

With one hand, press the button on the armrest and, with the other, lift the armrest. Then press the armrest down to the appropriate position to secure it. (images 10 and 11)


Image 10


Image 11

## SWITCH BETWEEN MANUAL MODE AND ELECTRIC MODE

**Manual mode:** Turn the lever in the direction of the unlock signal on the left and right motors (image 12).

**Electric mode:** Turn the lever in the direction of the lock signal on the left and right motors (image 13).

When switching the wheelchair from manual mode to electric mode, ensure that the left and right switches are in the same mode. Otherwise, it may tip over.

Image 12

Image 13

## 4.2. Checks prior to using the wheelchair

Before driving, carry out the following daily checks:

- **Battery check:** Before using your chair for the first time, charge the battery for 8 hours. Are the batteries sufficiently charged? The green lights on the battery indicator should be on.
- **Control pad check:** With the control system turned off, check that the joystick is not bent or damaged, and that it returns to the centre when you push and release it.
- **Seat check:** Check that all cushions are in place.
- **Visually inspect** the wheelchair to ensure that the armrests, etc., are correctly positioned on the wheelchair and that all clamps are properly tightened.
- **Checking for trapped clothing:** When starting the wheelchair, check that clothing does not interfere with its operation (i.e., that clothing is not too long). Before using the chair, always check that your clothing or accessories do not come into contact with the wheels or any other moving and/or rotating parts in which they could become entangled.
- **Checking weather conditions:** In winter, batteries have a reduced capacity. During a period of light frost, the capacity of the batteries is approximately 75% of normal capacity. At temperatures below -5°C, this will be approximately 50%. This factor reduces the range.
- If the foot platform rubs against the ground during use, it could break and cause injury. Therefore, handle with care and replace the platform if it is damaged.

### 4.3 Wheelchair handling

The wheelchair user is fully responsible at all times for complying with the relevant local safety regulations and guidelines.

Electric wheelchairs are powered by a control pad.

1. Switch on the control unit.
2. Set the maximum speed limit.
3. Move the joystick in the direction you wish to travel.
4. Moving the joystick further forward will make the wheelchair travel faster.

#### **DANGER!**

- If you have any doubts, do not risk crossing the road: wait until you are certain it is safe to cross.
- Always cross the road as quickly as possible; other vehicles may appear.
- Ensure that there are no objects in the path that could get caught in the mechanism of your wheelchair or in the rear wheel tyres. This could cause the wheelchair to stop abruptly.
- Driving over drains or grates could cause the front or rear wheels of the wheelchair to get stuck, resulting in the wheelchair stopping abruptly.

#### **ADVERSE CONDITIONS:**

Please remember that when operating the wheelchair in adverse conditions, such as on slippery surfaces, it may experience reduced grip and traction.

#### **DANGER!**

We recommend that you take extra precautions in these conditions, especially on slopes, as the wheelchair may become unstable or slip, causing injury.

### 4.4 Curves

#### **DANGER!**

Do not attempt to turn at high speeds. If you need to turn sharply, reduce your speed using the joystick or the speed adjustment on the control panel. This is especially important if you are crossing a slope or travelling downhill. Failure to follow this procedure could cause the wheelchair to tip over.

### 4.5 Brakes and emergency stop

There are three ways to brake the wheelchair.

- The simplest and safest way to stop the chair is to release the joystick. This will cause the chair to come to a controlled stop.
- Pulling the joystick back will cause the chair to stop abruptly with rapid braking.
- You can also stop the chair by turning off the control system while the chair is still moving (not recommended).

#### **WARNING!**

The method of switching the control system off and on again should only be used in an emergency situation, as the braking action is very abrupt.

## 4.6 Obstacles and kerbs

- Never descend a kerb backwards.
- Do not attempt to climb or descend flights of steps or use escalators. This would be very dangerous. You could suffer serious injury and the wheelchair could be damaged. This wheelchair is designed to climb a single step or kerb.
- We recommend that users with upper body instability use restraint and support systems to keep the upper body in position when going up or down ramps, kerbs or other obstacles.

### WHEN GOING UP OR DOWN A KERB:

- Approach the kerb or step head-on, at a 90° angle (figure 2).
- Drive forward slowly and progressively.
- Stop the chair as soon as the front wheels touch the kerb.
- Apply enough power to the motors to lift the front of the wheelchair so that it rests on the kerb, then apply a little more power so that the middle wheels climb the kerb smoothly.
- Keep the joystick in the forward position as much as possible.
- The maximum height of a step, obstacle or kerb is **25 mm**.

### LOWERING KERBS:

- Slowly and carefully move the chair forward towards the kerb until the front wheels are on the edge, again at a 90° angle to the kerb.
- Lower the kerb as slowly as possible with the drive wheels. Do not stop the chair when descending the kerb. You will feel safer if you can lean back, but if you cannot, do not worry, the wheelchair is stable.
- As long as you remain seated in it, you will be safe.
- All electric seat options should be in the start position. You may need to adjust the electric leg rests to allow enough clearance when going up or down a kerb.
- We recommend using a lap belt to increase your sense of security when going down the kerb.


Figure 2

## 4.7 Manual mode of the wheelchair

To use the wheelchair in manual mode (to push the wheelchair), the motors must be disconnected using the clutch lever. This feature was developed for use exclusively by the companion, never by the user.

**TO USE THE WHEELCHAIR IN MANUAL MODE:**

Turn the handle in the direction of the unlock signal on the left and right motors.

Manual Mode

Electric Mode

**⚠ DANGER!**

- Never leave the user unattended if the chair is in manual mode.
- Never set the lever to manual mode when travelling on a slope. When the release lever is in the "manual" position, the automatic parking brake is deactivated.
- The automatic parking brake only works when the lever is in electric mode.
- When you no longer need to push the chair, you must immediately set the release lever to electric mode.
- Ensure you have complete control of the chair when activating manual mode.
- Before activating manual mode, ensure the chair is on level ground.
- Push the chair by the handles. These provide secure points for holding the chair from behind and preventing falls or tip-overs.

## 5. Battery and charger

### 5.1 Charger requirements

The charger is used to charge the battery. Do not use the electric wheelchair while it is charging.

#### TECHNICAL DATA OF THE CHARGER

**Input voltage:** 100-240 V

**Output voltage:** 24 V DC

**Output current:** 3.0 A

The level of protection against water ingress is **IPX1**.

## 5.2 Using the charger

### 5.2.1 Type 1 (direct charging of the battery)

To charge the battery, connect the charger to the power supply and to the plug on the battery box.

#### FOLLOW THE INSTRUCTIONS BELOW TO COMPLETE THE UPLOAD PROCESS:

**Step 1:** Ensure that the charger slot is not blocked.

**Step 2:** Ensure that the electric wheelchair is switched off.

**Step 3:** Unplug the plugs connecting the battery box and the controller.

**Step 4:** Connect the charger's output plug to the battery box's power plug.

**Step 5:** Connect the charger's main plug to the power supply and the red light will turn on. A full charge takes between 8 and 10 hours; do not overcharge for more than 24 hours.

**Step 6:** To avoid shortening the battery's life, charge it at least once a month when the wheelchair is not in use.

### 5.2.2 Type 2 (charging via controller)

To charge the battery, connect the charger to the power supply and to the controller socket.

#### FOLLOW THE INSTRUCTIONS BELOW TO COMPLETE THE UPLOAD PROCESS:

**Step 1:** Ensure that the charger slot is not blocked.

**Step 2:** Ensure that the electric wheelchair is switched off.

**Step 3:** Ensure that the battery and controller are connected.

**Step 4:** Connect the charger's output plug to the socket located below the controller.

**Step 5:** Plug the charger's main plug into the power supply and the red light will turn on. A full charge takes between 8 and 10 hours; do not overcharge for more than 24 hours.

**Step 6:** To prevent shortening the battery's life, charge it at least once a month when the wheelchair is not in use.

### 5.2.3 Type 3 (indirect battery charging)

To charge the battery, you can remove it from the chair and take it wherever you wish to charge it.

#### FOLLOW THE INSTRUCTIONS BELOW TO COMPLETE THE UPLOAD PROCESS:

**Step 1:** Ensure that the electric wheelchair is switched off.

**Step 2:** Remove the battery from the wheelchair. Press and hold the two orange buttons and pull outwards.

**Step 3:** Ensure that the charger slot is not blocked.

**Step 4:** Connect the power plug from the battery box to the power supply and the red light will turn on. A full charge takes between 8 and 10 hours; do not overcharge for more than 24 hours.

**Step 5:** To avoid shortening the battery life, charge it at least once a month when the wheelchair is not in use.

Do not interrupt charging until the process is complete. Repeated use of a battery that is not fully charged will shorten its service life, so the battery should be charged to its maximum capacity whenever possible.

When the battery is fully charged, the power indicator will turn green. Do not interrupt charging before it is fully charged.

Once charging is complete, switch off the power supply, otherwise the battery will slowly discharge. Do not charge for more than 24 hours. Overcharging is dangerous.

**Users must follow the following rules to avoid hazards during charging:**

The electric wheelchair does not include a charger; use a standard national charger with an output voltage of 24 V/DC 1.5 A~5 A.

There must be good ventilation during charging. Do not expose the wheelchair to sunlight or humid environments.

The charging environment temperature range is 10°C to 50°C. If the ambient temperature range is exceeded, the battery will not function properly and will be easily damaged.

It is normal for the fan to make noise during charging. This is to cool the charger, so there is no need to worry.

Prevent liquid from entering the charger during charging. Do not place the charger on flammable objects such as fuel, footrests or seat cushions.

Stay away from flames while charging the battery. Flames can cause a fire or explosion of the battery.

Charging generates hydrogen; do not smoke while charging.

Do not unplug the power supply if the plug or your hands are wet, as this could cause an electric shock.

In the event of an unexpected accident and the user is injured, do not use or sit in the electric wheelchair while it is charging.

### 5.3 Battery use and maintenance

Incorrect operation when replacing the battery may cause a risk of explosion. Only replace it with a battery of the same type or the recommended type. Ensure that the battery terminals are correct. Key points for extending battery life: charge it frequently to keep it fully charged. It is best to fully charge the battery if the wheelchair is not in use. If it is not used for a long period of time, it is best to charge it twice a month.

## 6. System diagnosis

When the LED indicators flash, it means that the wheelchair has malfunctioned. Malfunctions can occur in the following parts: motors, brakes, battery, cable connections, etc.

By consulting the product's internal information, the nature of the malfunction can be detected through the diagnostic signal. The malfunction can be detected without the need for other service tools.

**AUDIO SIGNAL INDICATION**

Description of LED light	The meaning of LED light	Explanation and solution
All LED lights are off and no sound can be heard.	The power supply is disconnected, the wheelchair is in standby or sleep mode. The power supply has a poor connection. The fuse has tripped or blown.	
All the LED lights are on.	The power supply is connected and the self-diagnosis is working; the electric wheelchair can function properly.	The fewer LED lights that are on, the less battery power remains.
The red LED light on the far left is lit.	The battery power is extremely insufficient.	Charge immediately. Otherwise, the battery will not function properly and cannot be charged.
Two short beeps, the LED lights flash twice.	The left motor is not functioning correctly.	The left motor is incorrectly connected or the cable is disconnected.
Three short beeps, the LED lights flash three times.	The left magnetic brake is not functioning correctly.	The left magnetic brake is poorly connected or the cable is disconnected.
Four short beeps, the LED lights flash four times.	The right engine is not functioning correctly.	The right motor is incorrectly connected or the cable is disconnected.
Five short beeps, the LED lights flash five times.	The right-hand magnetic brake is not functioning correctly.	The right magnetic brake is incorrectly connected or the cable is disconnected.
Six short beeps, the LED lights flash six times.	The controller is in overcurrent protection mode.	Check the brakes and see if the motor drive mechanism is jammed. Check the current with an ammeter; if it is not excessive, the controller may be malfunctioning.
Seven short beeps, the LED lights flash seven times.	The joystick is not functioning correctly.	The joystick does not reset or the connector is loose.
Eight short beeps, the LED lights flash eight times.	The controller is not functioning correctly.	Consult your dealer for maintenance.
Nine short beeps, the LED lights flash nine times.	The controller is not functioning correctly.	Consult the service centre for servicing.

## 7. Controller

### 7.1 Control panel

### 7.2 Use of the Controller

#### POWER SWITCH

Press this button and the battery indicator lights will turn on from left to right. Press it again and all LED lights will turn off.

In case of emergency, you can switch off the device directly by pressing the power button.

#### SLEEP MODE

If the joystick is not used for more than 20 minutes, it will automatically switch off and the system will enter sleep mode. To activate it, press the power button.

## SPEED ADJUSTMENT

Depending on user habits and circumstances, the driving speed of the wheelchair is adjustable. Adjust the speed by pressing the decrease button or the increase button.

The speed is divided into five sections, ranging from 20%, 40%, 60%, 80% and 100% of the maximum speed.

First speed gear: 20% of maximum speed.

Fifth speed gear: 100% of maximum speed.

## HORN BUTTON

Press the horn button and it will sound until you release it.

## BATTERY CHARGE INDICATOR

After switching on the device, the battery indicator starts working. The battery indicator also shows the remaining battery capacity.

As shown in the image, the battery is fully charged.

When only the red or yellow LED light is on, the battery needs to be charged.

And for a long-distance trip, the battery must be fully charged. If only the red LED light is on, the battery is very low and users should charge it as soon as possible.

## USING THE JOYSTICK

The direction of movement of the wheelchair is controlled by a joystick. The range of movement of the joystick also controls the speed of movement.

Switch the device on or off when the joystick is in the centre position. Otherwise, the controller will emit an incorrect operating signal. Release the joystick to the centre position and the error will disappear. If the error persists, this part may be faulty. Do not use it and contact your dealer for repair.

## 8. Maintenance

Wheelchair maintenance includes cleaning, checking the wheels and battery, and charging the battery. For more thorough maintenance, contact your dealer.

### 8.1 Battery cleaning and maintenance

#### CLEANING

- Clean the wheelchair regularly.
- Clean parts that come into frequent contact with the user's body (such as the cushion, armrests, and controller) with a clean, slightly damp cloth. Dry them thoroughly.
- Avoid washing the wheelchair with water. The surface of the chassis has a protective coating that can be easily cleaned with a damp cloth.
- Do not use chemicals on the vinyl seat and armrests, as they may cause slipping or cracking. Use a damp cloth with neutral soapy water and dry them thoroughly.
- Also, do not use greasy products that may degrade the materials or make them slippery.
- When not in use, store the product in a dry place at room temperature.

#### WHEELS

Regularly check the condition of your tyres. When the tread depth is reduced to 1 mm, replace them with new ones.

#### BATTERY

Ensure that the battery is fully charged frequently. To prolong the battery's lifespan, we strongly advise users not to charge it until it has been completely depleted. If the battery indicator shows red (low battery), charge it immediately.

#### AFTER USE:

- Disconnect the power supply.
- Do not allow children or unconscious persons to use the wheelchair.
- Store the wheelchair at normal temperature to maintain its performance.

## 8.2 Wheelchair faults and checks

If a fault occurs during wheelchair operation, switch it off before checking it.  
Symptom: total loss of power and all LED lights on the controller panel off.

### VERIFICATION STEPS:

**Step 1:** Check whether the controller plug is loose.

**Step 2:** Check whether the connection between the controller plug and the battery box is loose. Reinsert the plug connector (hold the plug when removing it. Do not pull on the cable to avoid unnecessary damage to the cable). If, after performing the above checks, the wheelchair still does not regain power, or if users have any questions about the above checks, please contact your distributor.

The controller has a diagnostic system for monitoring the controller and motor. The controller indicates any malfunction of these parts. For further details, refer to the chapter on audio signal indication.

## 8.3 Maintenance check

The following is a checklist. The electric wheelchair should be checked in order according to our suggestions. When getting in or out of the wheelchair, some automatic checks are performed automatically. For your attention, we have specifically listed these self-check items in area A.

<b>Area A</b>	Before use, check that the following parts are correct: 1. Backrest. 2. Armrests. 3. Controller position. 4. Footrests. 5. Battery power. 6. Clutch/adjustment lever for switching between manual and electric mode.
<b>Area B</b>	Check the following parts monthly to prevent original parts from being lost or worn out: 1. Screws. 2. Brakes. 3. Clutch/adjustment lever for switching between manual and electric mode. 4. Front and rear wheels and tread depth. 5. Controller and charger connectors.
<b>Area C</b>	For safety reasons, it is necessary to carry out general maintenance every six months.

## 9. Warranty

The warranty is valid for three years (six months for batteries) from the date of the sales invoice, provided that the product has not been altered from its original configuration.

Within the warranty period, we will remedy any manufacturing defect at no cost, either by repair, replacement of parts or supply of a new product, at our discretion. To opt for repair or replacement, the owner of the product must go to the retailer, who will be solely responsible for handling the warranty with the manufacturer.

This warranty shall not be valid in cases of improper use or handling of the product, deterioration of the product due to external agents such as water or other harmful chemicals, obstructive or corrosive products, due to improper maintenance, lack of cleaning, loss of components or the use of non-original spare parts from our brand. For more information, please refer to our general terms and conditions of sale.

### 9.1 Warranty issues

#### WARRANTY CONTENTS

We have carefully designed this wheelchair for you. If there are any materials or manufacturing defects, we provide free repairs and lifetime maintenance in accordance with the time and conditions of the warranty cards.

#### THE WARRANTY DOES NOT COVER THE FOLLOWING CONDITIONS

Subjective impression, no operational problems.  
 Problems caused by lack of care or maintenance, or inadequate or incorrect maintenance.  
 Problems caused by impact or misuse of the product.  
 Wear and tear due to use and ageing (coating and plating surface, discolouration, etc.).  
 Operation different from that indicated in our manual or overload.  
 Unauthorised modifications.  
 External factors such as: soot, pharmaceutical products, bird droppings, acid rain, flying stones, metal dust, etc.  
 Natural disasters such as: typhoons, floods, fires, earthquakes, etc.

#### THE FOLLOWING FEES ARE NOT INCLUDED IN THE PRODUCT PRICE

Replacement of consumable materials, such as wheels, tyres, battery, plastic parts, motors, lubricants, etc.  
 Fees for inspection, adjustment, cleaning, etc.  
 Fees for maintenance checks.  
 Unauthorised modifications.

#### ACCEPTABLE WARRANTY

If you require the warranty, bring the electric wheelchair, take photos of the labels, and present the purchase invoice. If the user cannot provide these documents, the service will not be provided.

#### PROTECTING THE ENVIRONMENT

In order to protect the environment, when disposing of the packaging, product or any of its parts, strictly comply with the legal regulations of your community. Do not dispose of them indiscriminately.

## 10. Warning

- If you notice any defects, anomalies or malfunctions, notify the retailer or distributor immediately so that they can replace or repair the product.
- If you are under the influence of any medication that may affect your ability to drive, you should not use the wheelchair.
- To use the wheelchair safely, the user must have adequate vision.
- The wheelchair cannot accommodate more than one person at a time.
- Do not allow children to ride in the wheelchair without supervision.
- Any serious incident related to the product must be reported to Sucesores de Manuel García, S.L. and to the competent authority of the Member State in which the user and/or patient are established.

## 11. Others

### CONDITIONS FOR TRANSPORT AND STORAGE

During transport and storage, the electric wheelchair must be positioned correctly in accordance with the label instructions.

- Transportation should avoid moisture and sunlight, and be kept away from heat sources.
- In case electrical parts are damaged due to moisture, avoid storing the wheelchair in the rain, outdoors, or in damp places.
- Storage conditions:
  - Ambient temperature: -40 °C~+55 °C;
  - Relative humidity  $\leq$ 80 %;
  - Atmospheric pressure: 86 kPa~106 kPa.

## Annex to the report: Manufacturer's declaration on EUT

### ANNEX TO THE REPORT: MANUFACTURER'S DECLARATION ON EUT

1	<b>Manufacturer's guide and declaration: electromagnetic emissions</b>		
2	<p>The ELECTRIC WHEELCHAIR is designed for use in the electromagnetic environment specified below.</p> <p>The customer or user of the ELECTRIC WHEELCHAIR must ensure that it is used in such an environment.</p>		
3	Emission test	Compliance	Electromagnetic environment: Guide
4	RF Emissions CISPR 11	Group 1	The ELECTRIC WHEELCHAIR uses radio frequency energy solely for its internal operation. Therefore, its radio frequency emissions are very low and are not likely to cause interference with nearby electronic equipment.
5	RF Emissions CISPR 11	Class B	The ELECTRIC WHEELCHAIR is suitable for use in all establishments, including domestic ones and those directly connected to the low-voltage public electricity supply network that supplies domestic buildings.
6	Harmonic emissions IEC 61000-3-2	Class A	
7	Voltage fluctuations / flicker emissions IEC 61000-3-3	Compliance	

<b>Recommended separation distances between portable and mobile RF communications equipment and the ELECTRIC WHEELCHAIR</b>			
<p>The ELECTRIC WHEELCHAIR is designed for use in an electromagnetic environment in which radiated radio frequency interference is controlled. The customer or user of the ELECTRIC WHEELCHAIR can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radio frequency communications equipment (transmitters) and the ELECTRIC WHEELCHAIR, as recommended below, in accordance with the maximum output power of the communications equipment.</p>			
Maximum nominal power of the transmitter W	Separation distance according to transmitter frequency m		
	150 kHz to 80 MHz $d = \left[\frac{3.5}{V_1}\right]\sqrt{P}$	26 MHz to 800 MHz $d = \left[\frac{3.5}{E_1}\right]\sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.018	0.0345
0.1	0.38	0.057	0.1095
1	1.2	0.18	0.345
10	3.8	0.57	1.095
100	12	1.8	3.45
<p>For transmitters with a maximum output power not included in the above list, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the transmitter frequency, where P is the maximum output power of the transmitter in watts (W) as specified by the transmitter manufacturer.</p> <p><b>NOTE 1</b> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.  <b>NOTE 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p>			

### Manufacturer's guide and declaration: electromagnetic immunity

The ELECTRIC WHEELCHAIR is designed for use in the electromagnetic environment specified below. The customer or user of the ELECTRIC WHEELCHAIR must ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Level of compliance	Electromagnetic environment: Guide
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be made of wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Transient / surge electrostatic IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	The quality of the power supply should be that which is typical in a commercial or hospital environment.
Wave IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	The quality of the electrical network should be that of a typical commercial or hospital environment.
Voltage dips, interruptions and voltage fluctuations on power supply input lines IEC 61000-4-11	< 5% UT (>95% immersion in UT) for 0.5 cycles 40% UT (60% immersion in UT) for 5 cycles 70% UT (30% immersion in UT) for 25 cycles < 5% UT (>95% immersion in UT) for 5 seconds	< 5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds	The quality of the electrical network must be that which is normal in a typical commercial or hospital environment. If the user of the ELECTRIC WHEELCHAIR needs it to continue operating during power outages, it is recommended that it be powered by an uninterruptible power supply or a battery.
Mains frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	

**NOTE:** UT is the alternating current mains voltage before the test level is applied.

Manufacturer's guide and declaration: electromagnetic immunity			
The ELECTRIC WHEELCHAIR is designed for use in the electromagnetic environment specified below. The customer or user of the ELECTRIC WHEELCHAIR must ensure that it is used in such an environment.			
Immunity test	Test level IEC 60601	Level of compliance	Electromagnetic environment: Guide
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3V	<p>Portable and mobile RF communications equipment should not be used within the recommended separation distance, calculated from the equation applicable to the transmitter frequency, from any part of the ELECTRIC WHEELCHAIR, including cables.</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz a } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz a } 2.5 \text{ GHz}$ <p>where P is the maximum rated output power of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).<sup>b</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, must be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	20 V/m 26 MHz to 2.5 GHz	20V/m	
<p><b>NOTE 1</b> At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p><b>NOTE 2</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.</p>			
<p><sup>a</sup> The field strengths from fixed transmitters, such as mobile phone (cellular/wireless) base stations and land mobile radios, amateur radios, AM and FM radio broadcasts, and television broadcasts, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, consideration should be given to conducting an electromagnetic site survey. If the field strength measured at the location where the ELECTRIC WHEELCHAIR is used exceeds the applicable RF compliance level indicated above, the ELECTRIC WHEELCHAIR should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the ELECTRIC WHEELCHAIR.</p> <p><sup>b</sup> In the frequency range from 150 kHz to 80 MHz, field strengths should be below 3 V/m..</p>			

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Electric wheelchair  
Model: 9004 G-MOVE C2  
Instruction manual version 2.0  
Edition 04-2026

Jerry Medical Instrument (Shanghai) Co., Ltd.  
Room 116, Block A, Building 1, No. 288 Jiachang Road,  
Jiading District, Shanghai, China 201800

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MedPath GmbH