



PLATEADO
SILVER

MEDIDAS
MEASURE



1431SR	41 - 51	116	72 - 82	122	49	53	30	46 - 51	125	40
1431SRA			60 - 70			58	50			42
cm.									Kg.	

POSITRON

1431SR

SILLA DE POSICIONAMIENTO BASCULANTE Y RECLINABLE
TILTING AND RECLINING REHAB WHEELCHAIR



HASTA 125KG.
UP TO 125KG.



ALUMINIO
ALUMINIUM

1431SR

1431SRA

Estructura Frame	Reclinable de 90° a 135°. Basculante de 0° a 45°. Ángulo del asiento a 4° y ergonómico. Altura ajustable. Mecanismo de reducción de efecto cizalla durante la reclinación. Tilting from 90° to 135°. Reclining from 0° to 45°. Seat angle at 4°. Ergonomic seat. Adjustable height. Shear effect reduction mechanism during reclining.	
Ruedas Delanteras Front Wheels	180mm macizas. 180mm PU.	
Ruedas Traseras Rear Wheels	600mm macizas. 600mm PU.	300mm macizas. 300mm PU.
Reposabrazos Armrests	Ajustable en alto y ancho. Extraíble y abatible. Adjustable in height and width. Folding and detachable.	
Reposapiés Footrests	Elevable. Regulable en altura. Abatible interior y exterior. Elevating footrest, height adjustable. Swing in and out.	
Respaldo Backrest	Acolchado y anatómico. Reposacabezas estándar incluido. Padded and anatomically shaped. Headrest included.	
Accesorios Accessories	Taco abductor, soportes laterales. Leg abductor, trunk lateral supports.	



AZUL OSCURO
DARK BLUE

MEDIDAS
MEASURE



1432SR	41	123	63	31	126	50	80	40	120	29,5
	45		68							
cm.									Kg.	

BASIC II

1432SR

SILLA RESPALDO RECLINABLE
TILTING WHEELCHAIR



RECLINABLE
TILTING



HASTA 115KG.
UP TO 115KG.



ACERO
STEEL

Estructura Frame	Acero, plegable. Reclinable hasta 180°. Steel. Tilting up to 180°.	
Ruedas Delanteras Front Wheels	200mm macizas. 200mm PU.	
Ruedas Traseras Rear Wheels	600mm macizas MAG. 600mm PU MAG.	
Reposabrazos Armrests	Abatible y extraíble. Folding and detachable.	
Reposapiés Footrests	Elevables, ajustables en altura. Abatibles hacia el exterior. Elevating footrests, height adjustable. Swing away.	
Respaldo Backrest	Respaldo largo con cabezal incluido. Large backrest with headrest included.	
Accesorios Accessories	Anti-vuelcos. Reposapiés elevables. Porta oxígeno. Porta suero. Reposacabezas. Porta bastones. Anti-tipping wheels. Elevating footrests. Oxygen holder. Serum holder. Headrest. Cane holder.	